



GROWON MEDICAL TECHNOLOGY

DIVISION OF ELECTRO-CARE Since 1987

One Stop Solution

SURGICAL LASER



CORA 5 00



CORA SEAL



Cora Series

CORA-300



Light weight Compact easy to use, low cost solid state analog electro surgical unit with 250 W cut, coagulation Monopolar & bipolar mode.

MODE	DAIL	POWER	LOAD	OPERATING FREQ.
CUT	10	100 W	500	480 Khz
COAG	10	80 W	300	480 Khz
BIPOLAR	10	250	500	480 Khz

STANDARD ACCESSORIES



Chunk Handle with Electrodes Autoclavable



Bipolar Forcep with Cable Autoclavable



Monopolar Double Pedal Foot switch / Bipolar Single Pedal Foot switch



Disposable Hand Switch Pencil - Imported



Patient Plate with Cord Autoclavable

Workstation



OT LIGHT



Syringe pump



OT TABLE



ECG



ICU Monitor



ICU Ventilator



Z TROLLEY



FOGGER



SUCTION


DIVISION OF ELECTROCARE Since 1987
DIVISION OF ELECTROCARE Since 1987

MANY MORE EQUIPMENTS

DIVISION OF ELECTROCARE Since 1987



NICU Warmer



LED PHYSIOTHERAPY

Bubble CPAP



Hemodialysis Machine



MOBILUX-5



SS-2000



SHOCKWAVE PRO-1000



LASER CLASS IV



C Arm



BISTOS Fetal Monitor



DEFIBRILLATOR



GROWON MEDICAL TECHNOLOGY

- Surgical Laser Cautery & Vessel Sealing Systems
- ICU & Anesthesia Ventilator ■ Physiotherapy Equipments

 **7067785550,8817755102**

 E-mail : sales@grownmedical.com

 Website : www.grownmedical.com

 28-G, Industrial Area, Bagroda, Bhopal

 Branch Office: Plot No. 26, Sector-1 Shakti Nagar, Bhopal (M.P.)-24

