




Bili Flux NEO 200 - Series

*Treat Neonatal Jaundice with Confidence*

- LED Phototherapy with 18 high bright super LEDs, LED lamps specially made for Jaundice treatment.
- LED life time is **one lakh** burning hours i.e, more than 25 times compared to any other conventional phototherapy lamps.
- Treatment irradiance at skin level up to  $60\mu\text{w}/\text{cm}^2/\text{nm}$  at 30cm with wave length 420 to 500nm.
- Digital timer for total lamp usage and patient exposure.
- Only 18 watts power consumption, up to 80% power saving.
- Up to 30% faster serum bilirubin breakdown compared to conventional phototherapy.
- Optical design ensures uniform light distribution to the patient surface area.



**LED Phototherapy Stand  
Model: NEO 200**