

*Proven Excellent
Clinical Performance*

Modes of Operation

SERVO SKIN AIR PRE-WARM MANUAL

- > Micro processor based servo controller
- > Skin, Pre-warm & Manual control (Air mode optional)
- > Two bright digital displays combination of RED and GREEN continuous monitoring of baby temp and set temp, set temp easily convertible from °C to °F
- > Baby tray with high standard mattress, four sides clear acrylic walls, with front and both the sides collapsible, baby tray tilting facility both head up & head down, clear acrylic base for under surface Phototherapy.
- > German make ceramic heater with parabolic reflector to maintain uniform temp across the mattress with swiveling source box easy to take X-ray and procedures



Radiant Warmer Stand
Model : NEO 320



Radiant Warmer with Fixed Baby Cradle and Lower Drawers
Model : NEO 310



Infant Trolley
Model : NEO 600



Radiant Warmer with Fixed Baby Cradle
Model : NEO 300